

## MULTI-STIR INCUBATION BATH MARK III

D-IBMKIII

09/03

### EQUIPMENT SUPPLIED:

1. Water bath with rim to hold dispenser bottle and shaped to accommodate an IKA Labortechnik KMO 2 basic magnetic stirrer and to hold a polypropylene test-tube holder (as supplied).
2. One polypropylene test-tube rack to fit into water bath.
3. Sixteen 15mm magnetic stirrer bars.
4. Twenty, 20 x150 mm rimless glass test-tubes to fit supplied test tube rack.
5. Two pieces of corflute plastic to act as insulators between water bath and top of magnetic stirrer.

### EQUIPMENT NOT SUPPLIED :

1. IKA Labortechnik KMO 2 basic magnetic stirrer (use at setting 350).
2. Adjustable volume dispenser fitted to a one litre bottle.
3. Precision immersion heater similar to a Julabo PC.
4. IKA YellowLab Test Tube Shaker TTS. We strongly recommend this particular test-tube stirrer.

